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
## ENFit™ Lopez Valve™ For use with Gastric or Feeding Tubes

Non-Sterile. Use Aseptic Technique. Single Use Only. Device does not contain natural rubber latex. Device does not contain DEHP.

### DESCRIPTION

The ENFit Lopez Valve is a disposable, single use, three-way valve for use in gastric irrigation and aspiration. The valve contains a medication port to attach an enteral syringe that allows the administration of patient medication without breaking the line connection.

### INTENDED USE / INDICATIONS FOR USE

The ENFit Lopez Valve is intended for gastric irrigation and aspiration. The ENFit Lopez Valve is not intended for venous related use.

### TARGET PATIENT POPULATION / USER GROUPS

The ENFit Lopez Valve is intended to be used and handled by trained clinical personnel only. Clinical settings may include but are not limited to out-patient infusion centers, in-patient units, or emergency departments.

The patient population are those hospitalized patients in need of gastric irrigation or aspiration. The range of diagnosis is broad.

### CONTRAINDICATIONS

None known.

### WARNINGS – RESIDUAL RISK

Inspect device and read Instructions for Use prior to usage. A damaged device or improper use may lead to the following hazardous situations: delay of therapy, ergonomic stress, embolism from particulate, exposure to allergenic substances/biological contaminants/toxins/caustic substances. Harms potentially associated with these hazardous situations include reversible injury, permanent injury, or death.

### PRECAUTIONS

- Device is intended to be used and handled by trained clinical personnel only.
- Device is intended to connect to ISO 80369 compatible luer connectors.
- Do not use on intravenous lines. Device is compatible with enteral tubing luers and enteral syringes. Handle device in a similar nature as other single use enteral products.
- Device is not intended for use more than 72 hours.
- Inspect the device prior to each use for signs of wear (i.e. blocked luer, black lettering worn off, cracks). Replace if visible wear is observed.
- Tighten and secure all non-bonded connections and inspect for leaks before use.
- Device may be cleaned using a mild detergent and water when appropriate.
- Do not Reuse. Reuse negatively impacts performance potentially resulting in product failure/ contamination. Reuse may result in infections and/or contamination of tube feedings, which could result in gastrointestinal problems, such as diarrhea, bloating, etc.
- Change device in accordance with accepted standards in enteral therapy.
- Dispose according to Federal, State, and local regulations for treated infectious waste material.

### INSTRUCTIONS FOR USE

#### For Feeding:

1. Securely attach the male luer (A) on ENFit Lopez Valve to pump tubing. Tighten by hand only.
2. Attach female luer (B) to patients naso-gastric (NG), gastrostomy (G), or jejunostomy (J) tube. Tighten by hand only.
3. Attach pump set connector to male luer (A). Tighten by hand only. Ensure turn handle is OFF to the medication port (C) and end cap (E) is secure. Infuse as required.

#### For Suction:

1. Attach suction tubing to male luer (A). Attach female luer (B) to patient's naso-gastric (NG) tube.  
**Note:** adapter is required
2. Ensure turn handle is OFF to the medication male luer (C). Aspirate as required.

#### For Medication Delivery:

1. Attach syringe tip firmly into medication male luer (C). Tighten by hand only.
2. Turn stopcock handle until the "OFF" position points to male luer (A). Infuse or aspirate as required.
3. Return to suction or feeding by turning the stopcock handle until the "OFF" position points to the medication male luer (C).
4. Flush valve per facility protocol following administration.

### REMOVAL

1. Ensure return of stopcock handle in "OFF" position to the medication male luer (C).
2. Cap the medication male luer (C) when not in use.
3. Discard in accordance with facility protocol.

### CLINICAL BENEFITS

The ENFit Lopez Valve utilize current state-of-the-art technology. The featured benefits include:

- Maintains a closed system throughout enteral feeding, drug administration, and suction procedures.

- Provides access to enteral systems without opening lines, protecting from exposure to gastric fluids and bloodborne pathogens.
- Prevents inadvertent connection to vascular, respiratory, or epidural lines.
- Facilitates gastric lavage without risk of splashing
- Design conforms to ISO 80369 standard and guidelines for enteral feeding connectors

Use of the ENFit Lopez Valve is consistent with the current knowledge and state of art for devices used in gastric or enteral therapy. In normal use conditions, the clinical benefits of ENFit Lopez Valve outweigh the associated risks.

### ADVERSE EVENTS

None known

### INCIDENT REPORTING

Serious incidents that occur during use of the product should be reported to ICU Medical. These reports can be submitted at [www.icumed.com](http://www.icumed.com). Incidents should also be reported to the relevant local regulatory authority where required.

These products are covered by one or more U.S. and foreign patents.

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