Lopez Valve™


ICU Medical B.V.

 Hofspoor 3, 3994 VZ Houten
 The Netherlands

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Lopez Valve™ Closed Enteral Tube Valve

Non-Sterile. Use Aseptic Technique. Single Use Only. Device does not contain natural rubber latex. Device does not contain DEHP.

DESCRIPTION

The Lopez Valve is a disposable, single use, three-way valve for use in gastric irrigation and aspiration. The valve contains a medication port to attach an enteral syringe that allows the administration of patient medication without breaking the line connection. The Lopez Valve is not intended for venous related use.

INTENDED USE / INDICATIONS FOR USE

The Lopez Valve is intended for gastric irrigation and aspiration. The medication port of the Lopez Valve allows the healthcare worker to give the patient medication without breaking the line connection. The Lopez Valve provides a closed system that prevents the exposure of the healthcare worker to gastric fluids or blood-borne pathogens.

TARGET PATIENT POPULATION / USER GROUPS

The Lopez Valve is intended to be used and handled by trained clinical personnel only. Clinical settings may include but are not limited to out-patient infusion centers, in-patient units, or emergency departments.

The patient population are those hospitalized patients in need of gastric irrigation or aspiration. The range of diagnosis is broad.

CONTRAINDICATIONS

None known.

WARNINGS – RESIDUAL RISK

Inspect device and read Instructions for Use prior to usage. A damaged device or improper use may lead to the following hazardous situations: delay of therapy, ergonomic stress, embolism from particulate, exposure to allergenic substances/biological contaminants/toxins/caustic substances. Harms potentially associated with these hazardous situations include reversible injury, permanent injury, or death.

PRECAUTIONS

- Device is intended to be used and handled by trained clinical personnel only.
- Do not use on intravenous lines. Device is compatible with enteral tubing luer and enteral syringes. Handle device in a similar nature as other single use enteral products.
- Device is not intended for use more than 72 hours.
- Inspect the device prior to each use for signs of wear (i.e. blocked luer, black lettering worn off, cracks). Replace if visible wear is observed.
- Tighten and secure all non-bonded connections and inspect for leaks before use.
- Do not Reuse. Reuse negatively impacts performance potentially resulting in product failure/ contamination.
- Change Lopez Valve in accordance with accepted standards in enteral therapy.
- Dispose according to Federal, State, and local regulations for treated infectious waste material.

INSTRUCTIONS FOR USE

For Feeding:

1. Securely attach the universal adapter (D) to short-tapered port on Lopez Valve (A) and insert pump tubing into adapter until secure. Tighten by hand only.
2. Attach longest tapered port on Lopez Valve (B) to patient's naso-gastric (NG), gastrostomy (G), or jejunostomy (J) tube. Tighten by hand only.
3. Ensure turn handle is OFF to the medication port (C) and end cap (E) is secure. Infuse as required.

For Suction:

1. Attach suction tubing to male luer (A). Attach female luer (B) to patient's naso-gastric (NG) tube.
Note: universal adapter (D) is not required
2. Ensure turn handle is OFF to the medication port (C) and end cap (E) is secure. Aspirate as required..

For Medication Delivery:

1. Insert appropriate syringe tip firmly into medication port on Lopez Valve (C), twist until secure.
2. Turn stopcock handle until the OFF position points to the short-tapered port (A). Infuse as required.
3. Flush valve per facility protocol following administration.

REMOVAL

1. Ensure return of stopcock handle in OFF position to the medication port (C).
2. Cap medication port (C) and suction/adaptor port (A) when not in use (F).
3. Discard in accordance with facility protocol.

CLINICAL BENEFITS

The Lopez Valve utilize current state-of-the-art technology. The featured benefits include:

- Maintains a closed system throughout enteral feeding, drug administration, and suction procedures.
- Eliminates the use of nasogastric (NG) tube plugs and poor syringe connections while keeping the clinician safe from accidental exposure to infectious bodily fluids.

- Facilitates gastric lavage without risk of splashing

Use of the Lopez Valve is consistent with the current knowledge and state of art for devices used in gastric or enteral therapy. In normal use conditions, the clinical benefits of Lopez Valve outweigh the associated risks.

ADVERSE EVENTS

None known

INCIDENT REPORTING

Serious incidents that occur during use of the product should be reported to ICU Medical. These reports can be submitted at www.icumed.com. Incidents should also be reported to the relevant local regulatory authority where required.

These products are covered by one or more U.S. and foreign patents.

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